



## Special Fireproof Mattress / Special Fireproof Pillow

Järven's special fireproof mattress and special fireproof pillow are for extremely demanding and fire hazardous environments. The products fulfil all fire standards in the world that are known to us. They are easy to clean and disinfect, as well as extremely durable and fluid proof.

The fire and smoke emission hazards are minimized by the unique cover which is resistant to open gas burner flames and glowing particles. The products contain a fabric that consists of Aramid, a material used in

bulletproof vests that can withstand extremely high temperatures. The coating made of aramid-based composite materials is high frequency welded, which makes the mattress and pillow exceptionally durable and liquid-proof.

### Warranty

1 year for manufacturing defects. NOTE: The warranty does not apply if damage to the mattress or pillow is caused by mechanical impact from pens or other sharp objects.



Please read the instruction manual before using the product.

### Fire protection

Fulfills all fire standards in the world that are known to us. The mattress core 37 kg / m<sup>3</sup> is made of polyether, which meets the fire requirements NT-Fire 055. The inner pillow is made of 100% Trevira CS, a flame retardant textile with achievable flammability class B1/M1. The cover is resistant to open gas burner flames and glowing particles.

Fire standards achieved:

SS 876 00 10 test method NT-Fire 055, BS 6807:2006 Crib 7, FAR 25.853

### Item no.

32560  
32570  
32500  
32510  
32520  
32530  
32540  
32550  
32580  
32929

### Size

190-205 x 70 x 10 cm  
190-205 x 75 x 10 cm  
190-205 x 80 x 10 cm  
190-205 x 85 x 10 cm  
190-205 x 90 x 10 cm  
190-205 x 105 x 10 cm  
190-205 x 110 x 10 cm  
190-205 x 120 x 10 cm  
190-205 x 150 x 10 cm  
Pillow 50 x 60 cm

### User weight

Up to 130 kilograms.

### Recycling

The mattress is fully recyclable and should be sorted for energy recovery. No dismantling is required.



**Järven Health Care**

Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN  
Telephone: +46 (0)660-22 28 50 | Telefax: +46 (0)660 - 824 10  
www.jarven.se | healthcare@jarven.se