



Järven CPR-board

Safety · Strenght · Hygiene



Järven's new and sturdy, hygienic CPR-board is manufactured in hard plastic and is developed in close cooperation with experienced CPR-personnel. The grip friendly handles and the soft bevelled edges give you an effective CPR-board when needed

Care instructions

Järven CPR-board is easy to clean and disinfect, reducing the risk for infections. For normal cleaning, it is recommended to clean with soap solution and water. When disinfecting you can use an alcohol based cleaning agent with surfactant.

Item no.

19-000000-001

Size

70 x 50 x 1 cm



Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN
 Phone: +46 660-22 28 50 | Fax: +46 660 - 824 10
www.jarven.se | healthcare@jarven.se