

JÄRVEN OMLÄGGNINGSKUDDEN

A cushion for facilitating bandaging and other care procedures.



Järven Omläggingskudden Rev 2019.03

Omläggingskudden from Järven is designed to simplify daily care work. The part of the body which requires bandaging or care is placed on the cushion to provide an ergonomic work setting.

In the case of bandaging situations, it is often common that two people must assist – one to hold the limb and one to carry out the bandaging. With this support cushion a single person can manage both tasks. The shape of the cushion allows it to be rolled easily during the procedure. The cushion's surface layer is **wipeable for maintaining hygiene**.

Material

Filler material: Cold foam
Surface layer: Polyurethane (PU)

Labelling

Date of manufacture, EAN, care instructions & contact details.

Warranty

1 year for manufacturing defects.

Recycling

The cushion is completely energy-recoverable and can be sorted as combustible waste. No disassembly is necessary when disposing of the cushion.

Product care

Detailed information is contained in the product's user guide.

Item number

13-110000-001

Size

H: 27 cm x W: 29.5 cm

Weight

670 g

Design protection: DM/083438



Read the user guide before using the product.

