



## Järven Rescue Mattress



The rescue mattress can be used for quick and safe evacuation of people who are unable to evacuate themselves in case of fire or other hazard.

With the Järven safety mattress, one single person can with some simple manoeuvres evacuate a person in distress and it provides additional security in those cases an evacuation process can include injured, disabled and/or physically handicapped people.

The rescue mattress is symmetrically designed with an underside that slides as well on wet as on dry surfaces. It is X-ray transparent and equipped with Velcro straps that provide a firm fixation. The underside is heat reflective.



Buy an additional sign and increase the visibility!

Item no.	Description	Size (LxBxH)
11-020000-001	Rescue mattress	190 x 48 x 8 cm
11-020100-001	Sign + fastener	30x7 cm

**Warranty**  
1 year for manufacturing defects

**X-ray**  
The rescue mattress is X-ray transparent

**Labelling**  
The rescue mattress is labelled with evacuation instructions, care instructions and the manufacturer's contact details.

**Recycling**  
The rescue mattress is fully energy recoverable and should be sorted as combustible waste. No dismantling is required.



Please read the instruction manual before using the product.



Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN  
Telephone: +46 (0)660-22 28 50 | Telefax: +46 (0)660 - 824 10  
www.jarven.se | healthcare@jarven.se