

Evacuation



Järven Rescue Sheet

With Järven rescue sheets placed between the mattress and the bed, the patient can quickly be evacuated in case of fire or other hazard. One single person can evacuate a care recipient with a simple manoeuvre.

The Järven rescue sheet has an extra long pull handle at each short side that enable the staff to perform the evacuation with straight backs, thus reducing the physical strain.

The fastening straps contain strong Velcro and are a clear orange safety colour.

With a rescue sheet that meets the requirement of Swedish standard

SS 876 00 39:2012 one can feel safe with the knowledge that it is durable and can be used for evacuations in both dry and wet circumstances.

Item no. For mattress width

9900	78 - 84 cm
9901	70 - 77 cm
9902	85 - 90 cm
9903	100 - 105 cm



SS 876 00 39:2012 Rescue sheet standard

The Järven rescue sheet meets the requirements of Swedish standard SS 876 00 39:2012, and this guarantees, among other things, that the handles can withstand a load of 1,000 N (100 kg), which is the approximate equivalent of a patient weight of 210 kg. The sliding ability is high on both wet and dry surfaces.

Warranty

1 year for manufacturing defects

X-ray

The rescue sheet is X-ray transparent

Labelling

The rescue sheet is labelled with evacuation instructions in 8 languages, GS1 number, care instructions and manufacturer's contact information

Flame retardant

Meets the requirements of SS 876 00 01

Recycling

The rescue sheet is fully recyclable and should be sorted for energy recovery. No dismantling is required.



Please read the instruction manual before using the product.



Järven Health Care

Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN
Telephone: +46 (0)660-22 28 50 | Telefax: +46 (0)660 - 824 10
www.jarven.se | healthcare@jarven.se